



Everything you need to know, before you go to the solarium

PROTECT YOUR EYES

Always wear UV protective glasses when sunbathing.

MEDICATIONS

If you are taking medication, we recommend that you consult a doctor before sunbathing. Certain medications can make your skin particularly sensitive to sunlight.

SUNSCREEN

Do not use outdoor sun products when using the tanning bed. Thoroughly remove all cosmetics from your skin before your tanning session. Many cosmetics contain sunscreen. This will affect your tan. Want the best tanning results possible? Ask the tanning bed consultant for a special product suitable for use with the Ergoline tanning beds.

SUNBATHING AT THE HIGHEST LEVEL

With Ergoline sunbeds you can enjoy all the benefits of the sun. Every skin type is different. To achieve tan, it is important that you follow the advice of our qualified professionals. In the buildup phase, it is important to use the tanning bed 2 to 3 times a week, with the number of minutes and the type of tanning bed tailored to your skin type. Do you want to keep your tan? Then make sure you have a tanning session every one to two weeks. It's important not to put too much stress on your skin. Always take a break of at least 48 hours between tanning sessions.

FIRST SUN SESSION(S)

Which tanning bed and which tanning session is right for you depends on your skin type. Ask the Ergoline sun consultant for a customized consultation.

PLEASE NOTE! Do you recognize yourself with any of the following characteristics? Then you should not use the tanning bed. Please inform the sun consultant.

- Anyone under 18 years of age.
- Anyone with skin type 1.
- Anyone with or large number of moles.
- Anyone who suffers or has suffered from skin cancer.
- Anyone who has suffered burns from outdoor sun exposure. - Anyone who has experienced irritation or reddening of the
- skin after initial sun exposure.





SKIN TYPE 1



SKIN TYPE 2



SKIN TYPE 3



SKIN TYPE 4

Natural skin color	very pale
Freckles	very often
Natural hair color	reddish to reddish blond
Eye color	blue, grey
Sunburn	always and painfully
Tanning	none
Erythema effective irradiation	200 jm³

pale
often
blonde to brown
blue, green, grey to brown
almost always painful
barely to moderately
250 jm³

light to light brown
rare
dark blond to brown
grey, brown
rare to moderate
gradulally
350 jm³

light brown, olive	
none	
dark brown	
brown to dark brown	
rare	
fast and deep	
450 jm³	





